

Mini Wave Infinity **WATERPROOF**
The first, true wireless underwater EMG system in the world.



c@meta




miniwave ∞

wireless EMG, accelerometers and data logger

TECHNICAL SPECIFICATIONS

- Memory equipped probes (7 grams).
- Unbeatable signal quality (no transmission of data).
- Suitable for water sports and rehabilitation.
- 2.000 Hz sampling rate at 16 bits.
- Up to 16 EMG channels with accelerometers.
- Asynchronous data logging for up to 8 hours underwater.
- Works as a regular Mini Wave Infinity system out of water.
- Synchronization possible with Remote control and TTL trigger.



COMETA SRL

via Giovanni Falcone, 43 - 20010 Bareggio, Milano, Italy

T. +39 02 91410582 • M. sales@cometasystems.com

www.cometasystems.com